

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property Organization
International Bureau




(43) International Publication Date
13 September 2001 (13.09.2001)

PCT

(10) International Publication Number
WO 01/67108 A2

- (51) International Patent Classification⁷: G01N 33/68 (74) Agent: LUCAS & CO; 135 Westhall Road, Warlingham, Surrey CR6 9HJ (GB).
- (21) International Application Number: PCT/EP01/02894
- (22) International Filing Date: 12 March 2001 (12.03.2001)
- (25) Filing Language: English
- (26) Publication Language: English
- (30) Priority Data:
0005683.8 10 March 2000 (10.03.2000) GB
0006064.0 14 March 2000 (14.03.2000) GB
- (71) Applicant (for all designated States except US): UNIVER-SITE DE GENEVE [CH/CH]; Rue General-Dufour 24, Case Postale, CH-1211 Geneva 4 (CH).
- (72) Inventors; and
- (75) Inventors/Applicants (for US only): HOCHSTRASSER, Denis, Francois [CH/CH]; Chemin de Savonniere, Col-longe-Bellerive, CH-1245 Geneva (CH). SANCHEZ, Jean-Charles [CH/CH]; Chemin Frank-Thomas 42, CH-1208 Geneva (CH). ZIMMERMANN, Catherine, Gabrielle [CH/CH]; Rue Maunoir 48, CH-1207 Geneva (CH). GUILLAUME, Elisabeth [FR/FR]; 8A, rue du Capitaine Charles Dupraz, F-74100 Annemasse (FR).
- (81) Designated States (national): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZW.
- (84) Designated States (regional): ARIPO patent (GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).
- Published:**
— without international search report and to be republished upon receipt of that report
- For two-letter codes and other abbreviations, refer to the "Guidance Notes on Codes and Abbreviations" appearing at the beginning of each regular issue of the PCT Gazette.

(54) Title: DIAGNOSTIC ASSAY FOR TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHIES



(57) Abstract: Heart and brain fatty acid binding proteins (H-FABP, B-FABP) are markers for TSEs, especially CJD. The invention provides a diagnostic assay for either of these markers, preferably by enzyme immunoassay using a specific antibody thereto. Since H-FABP is also a marker for acute myocardial infarction (AMI), to distinguish CJD from AMI requires an assay specific to AMI, e.g. using troponin-1 or CK-MB as a marker, also to be carried out.



WO 01/67108 A2

- 1 -

"DIAGNOSTIC ASSAY FOR TRANSMISSIBLE SPONGIFORM
ENCEPHALOPATHIES"

BACKGROUND OF THE INVENTION

Field of the invention

5 This invention is in the field of diagnostic assay
using a protein or an antibody thereto.

Description of the related art

 Transmissible spongiform encephalopathies (TSEs) are
neurodegenerative diseases of the central nervous system.
10 They can be transmitted, inherited or occur sporadically
and are observed in animals, e.g. as bovine spongiform
encephalopathy (BSE) in cattle or scrapie in sheep, as
well as in humans as Creutzfeldt-Jakob disease (CJD),
Gerstman Sträussler Scheinker syndrome, Fatal Familial
15 Insomnia or Kuru. They have a long incubation period,
leading to ataxia, dementia, psychiatric disturbances and
death. Neuropathological changes include vacuolar
degeneration of brain tissue, astrogliosis and amyloid
plaque formation. The diseases are difficult to diagnose
20 pre-mortem.

 The cerebrospinal fluid (CSF) of CJD patients
displays two additional polypeptide by two-dimensional
polyacrylamide gel electrophoresis [Harrington, M.G. New
England Journal of Medicine 315, 279 (1986), Hsich, G.,
25 Kenney, K., Gibbs, C.J., Lee, K.H. & Harrington, M. B.
New England Journal of Medicine 335, 924 (1996).] The
function of these 14-3-3 polypeptides remain unclear in
TSE. They can be used in a pre-mortem test for CJD
diagnostic evaluation, but have low specificity.

30 Monoclonal antibodies to the abnormal form of prion
protein are available and can be used in an enzyme-linked
immunoassay, as described in PCT Specifications WO
98/23962 and 98/32710 and Schmerr, M.J., the Beckman

- 2 -

Coulter Pace Setter Newsletter 3(2),1-4 (June 1999), but these procedures have not yet been fully developed.

Development of new non-invasive blood CJD and BSE markers would help clinicians to establish early
5 diagnosis.

SUMMARY OF THE INVENTION

It has now surprisingly been found that two fatty acid binding proteins (FABP), known as heart (H-FABP) and brain (B-FABP), are markers for TSEs. Thus, the
10 invention provides a diagnostic assay for a TSE or the possibility thereof in a sample of body fluid taken from a subject suspected of suffering from the TSE, which comprises determining the concentration of heart or brain fatty acid binding protein (H-FABP or B-FABP) in the
15 sample. The method is especially applicable to the diagnosis of CJD, especially new variant CJD, in human patients, and to BSE in ruminant animals such as cattle.

Conveniently the method is carried out using an antibody to H-FABP or B-FABP, whereby the extent of the
20 reaction between the antibody and the FABP in the sample is assayed and related to the concentration of FABP in the sample. The concentration thus determined is used to make or assist in making a diagnosis.

The present invention enables an assay of high
25 sensitivity, specificity and predictive accuracy for CJD to be carried out. "Sensitivity" is defined as the percentage of true positives given by the assay on samples taken from patients in whom clinical examination has confirmed CJD. "Specificity" means the percentage of
30 true negatives given by the assay on control samples, i.e. from patients in whom clinical examination has not revealed CJD. "Predictive accuracy" means the ratio of true positives to total positives (true + false) expressed as a percentage.

- 3 -

H-FABP is a known marker of acute myocardial infarction (AMI), see Ishii, J. et al., "Serum concentrations of myoglobin vs human heart-type cytoplasmic fatty-acid binding protein in early detection of acute myocardial infarction", Clinical Chemistry 1997;43 1372-1378. Therefore, in order to use an assay for H-FABP for the diagnosis of CJD in humans to better advantage, it is desirable to perform another kind of assay for AMI (one in which the marker is not a FABP) in order to eliminate from the diagnosis for CJD those patients who are positive in the AMI assay.

Thus, in a particular embodiment, the invention provides a method which comprises determining the concentration of H-FABP in a first assay, as defined above, whereby a positive result indicates either a CJD or acute myocardial infarction, and which further comprises carrying out a second diagnostic assay, for acute myocardial infarction (AMI) only, whereby a positive result in the H-FABP assay and a negative result in the assay for AMI indicates that the patient might be suffering from CJD. Assays using Troponin-I and Creatine Kinase-MB (CK-MB) as early biochemical markers of acute myocardial infarction (AMI) are well known and suitable for the above purpose.

A similar H-FABP and also a brain-specific fatty acid binding protein (B-FABP) have been found in the brain of mice, see Pu, L. et al., Molecular and Cellular Biochemistry 198, 69-78 (1999). Brain H-FABP (not to be confused with B-FABP) is believed to differ from heart H-FABP by a single amino acid substitution. However, B-FABP differs considerably. Sellner, P.A. et al., "Development role of fatty acid binding proteins in mouse brain" Dev. Brain Res. 89, 33-46 (1995), estimated the DNA homology at 69%, while A. Schreiber et al.,

- 4 -

"Recombinant human heart-type fatty acid binding protein as standard in immunochemical assays", Clin. Chem. Lab. Med. 36(5), 283-288 (1998), mention 64% amino acid sequence homology and that a monoclonal antibody to human H-FABP is cross-reactive with human B-FABP to the extent of only 1.7%.

Now that the present inventors have found that H-FABP is a marker for CJD, it is a very reasonable prediction that B-FABP will also be. Since B-FABP is specific to brain tissue and does not appear to react significantly with a monoclonal antibody to H-FABP, it will not give positives for AMI, making a separate assay for AMI unnecessary.

BRIEF DESCRIPTION OF THE DRAWINGS

The Figure is a graphic representation on the y-axis of H-FABP concentration represented by optical density measurement at 405 nm, as determined by the method of the invention, for (a) a control group having neither CJD nor AMI (b) a group having AMI and (c) a group having CJD.

DESCRIPTION OF PREFERRED EMBODIMENTS

For the method of assay, the sample can be taken from any convenient body fluid of the subject, but preferably plasma or serum (rather than whole blood). Cerebrospinal fluid (CSF) is another useful fluid, particularly when testing animals such as cattle.

The method is considered applicable to all types of TSE, including those referred to above, and to any human or animal suffering or suspected of suffering therefrom. Particularly, the invention is applicable to all types of CJD in humans, including new variant, sporadic and genetic (familial). Further, it is applicable to BSE in cattle and BSE-like disease in other animals, e.g. deer.

The marker, H-FABP or B-FABP, is preferably measured by an immunoassay, using a specific antibody to H-FABP

- 5 -

and measuring the extent of the antigen (H-FABP or B-FABP)/antibody interaction. For the diagnosis of human patients, the antibody is preferably anti-human H-FABP or B-FABP. Similarly, if the subject is an animal the antibody is preferably anti- to the H-FABP or B-FABP of the same animal variety, e.g. anti-bovine H-FABP or B-FABP if the patient is bovine. However, there is some cross reactivity of the antibodies between species, often enabling a heterologous antibody to be used: for example anti-rat/mouse H-FABP can be used to detect BSE in cattle. It may be a monoclonal antibody (conveniently mouse) or an engineered antibody. Preferably a mouse anti-human, anti-bovine etc. monoclonal antibody is used. Antibodies to H-FABP are known, e.g. 66E2 and 67D3 described by Roos, W. et al., "Monoclonal antibodies to human heart type fatty acid-binding protein", J. Immunol. Methods 183 149-153 (1995). Antibody 66E2 is commercially available. Also, the usual Köhler-Milstein method may be used to raise H-FABP or B-FABP antibodies. The source of protein for this purpose can be the naturally derived or recombinant DNA-prepared protein. Recombinant human H-FABP and B-FABP have been described by Schreiber, A. supra and Shimizu, F. et al., "Isolation and expression of a cDNA for human brain fatty acid binding protein (B-FABP)", Biochim Biophys. Acta 1354, 24-28 (1997), respectively. Less preferably, the antibody may be polyclonal.

Any known method of immunoassay may be used. A sandwich assay is preferred. In this method, a first antibody to the FABP is bound to the solid phase such as a well of a plastics microtitre plate, and incubated with the sample and with a labelled second antibody specific to the H-FABP or B-FABP to be detected. Alternatively, an antibody capture assay could be used here, the test

- 6 -

sample is allowed to bind to a solid phase, and the anti-FABP antibody is then added and allowed to bind. After washing away unbound material, the presence or amount of antibody bound to the solid phase is determined using a
5 labelled second antibody, anti- to the first.

In another embodiment, a competition assay could be performed between the sample and a labelled FABP or a peptide derived therefrom, these two antigens being in competition for a limited amount of anti-FABP antibody
10 bound to a solid support. The labelled FABP or peptide could be pre-incubated with the antibody on the solid phase, whereby the FABP in the sample displaces part of the FABP or peptide thereof bound to the antibody.

In yet another embodiment, the two antigens are
15 allowed to compete in a single co-incubation with the antibody. After removal of unbound antigen from the support by washing, the amount of label attached to the support is determined and the amount of protein in the sample is measured by reference to standard titration
20 curves established previously.

The label is preferably an enzyme. The substrate for the enzyme may be colour-forming, fluorescent or chemiluminescent.

It is highly preferable to use an amplified form of
25 assay, whereby an enhanced "signal" is produced from a relatively low level of protein to be detected. One particular form of amplified immunoassay is enhanced chemiluminescent (ECL) assay. Here, the antibody is preferably labelled with horseradish peroxidase, which
30 participates in a chemiluminescent reaction with luminol, a peroxide substrate and a compound which enhances the intensity and duration of the emitted light, typically 4-iodophenol or 4-hydroxycinnamic acid.

- 7 -

Another preferred form of amplified immunoassay is immuno-PCR. In this technique, the antibody is covalently linked to a molecule of arbitrary DNA comprising PCR primers, whereby the DNA with the antibody attached to it is amplified by the polymerase chain reaction. See Hendrickson, E.R. *et al.*, *Nucleic Acids Research* 23, 522-529 (1995) or Sano, T. *et al.*, in "Molecular Biology and Biotechnology" ed. Robert A. Meyers, VCH Publishers, Inc. (1995), pages 458 - 460.

5 The signal is read out as before.

In a particularly preferred procedure, an enzyme-linked immunosorbent assay (ELISA) was developed to detect H-FABP in serum. Since H-FABP is a marker for AMI as well, Troponin-I or CK-MB levels were assayed in order to exclude any heart damage. As described in the Example, these assays were assessed in serial plasma and CSF samples, from patients lacking AMI and CJD, patients with AMI, patients with dementia and patients with confirmed CJD through autopsy. The sensitivity, specificity and predictive accuracy for H-FABP in CJD above a suitable cut-off level were all 100%. Thus, H-FABP detection combined with the Troponin-I or CK-MB assay provides a useful serum marker of CJD diagnosis or brain damage.

15

20

The use of a rapid microparticle-enhanced turbidimetric immunoassay, developed for H-FABP in the case of AMI, Robers, M. *et al.*, "Development of a rapid microparticle-enhanced turbidimetric immunoassay for plasma fatty acid-binding protein, an early marker of acute myocardial infarction", *Clin. Chem.* 44, 1564-1567 (1998), should drastically decrease the time of the assay. Thus, the full automation in a widely used clinical chemistry analyser such as the "COBAS" MIRA Plus system from Hoffmann-La Roche or the "AxSYM" system from

25

30

- 8 -

Abbott laboratories should be possible and applied for routine clinical diagnosis of CJD.

The H-FABP or B-FABP can be measured by other means than immunoassay. For example, the sample can be
5 subjected to 1 or 2-DE gel electrophoresis and the amount of the FABP estimated by densitometric scanning of the gel or of a blot therefrom.

The assay of the invention can be used together with one or more other pre-mortem assays for the TSE,
10 including specifically those assays described above. Such combined procedures are particularly useful in diagnosing BSE in ruminant animals such as cattle.

The following Examples illustrate the invention.

EXAMPLE 1

15 Materials And Methods

Patients

The study population consisted of 3 age-and-gender matched control patients (Control group), 3 confirmed AMI patients (AMI group), 3 confirmed dementia patients
20 (dementia group) and 3 confirmed CJD patients (CJD group). The Control group included 2 men, mean age 66, range 46-86 years, and 1 woman, age 63 years. The AMI group included 2 men, mean age 65, range 40-90 years, and 1 woman, age 72 years. The dementia group included 2 men,
25 mean age 65, range 43-87 years, and 1 women, age 64 years. The CJD group included 2 men, mean age 68, range 62-74 years, and 1 woman, age 65. Blood and CSF samples were collected for each patient of the CJD. Blood samples were collected in dry heparin-containing tubes. After
30 centrifugation at 1500g for 15min at 4°C, tubes were stored at -20°C until analysis. Patients from the CJD group underwent serial clinical evaluations by neurologists in order to confirm CJD diagnosis. Patients from the AMI group were admitted to the hospital with a

- 9 -

confirmed AMI (Troponin-I concentration $>2\text{ng/ml}$). A clinical evaluation was performed on all the patients from the control group to exclude CJD and AMI.

Measurement of brain and heart H-FABP

5 H-FABP levels were measured in plasma by a sandwich ELISA. A 96-well polystyrene microplate (NUNC) was coated with 100 microlitres/well goat anti-human FABP, detecting all isoforms (Spectral Diagnosis HC, Ontario, USA), 20 micrograms/ml in carbonate buffer 0.1M pH 9.6, overnight
10 at 4°C . The plate was automatically washed with PBS (15mM Na_2PO_4 -120mM NaCl-2.7mM KCl pH 7.4, Sigma) on a BioRad NOVAPATH™ washer. Every washing step was performed with fresh PBS. Non-specific binding sites were blocked with 200 microlitres/well 2% casein in carbonate buffer for 2h
15 at 37°C . After the washing step, the samples were pipetted in duplicate at 100 microlitres/well. The plate was incubated 2h at 37°C . After the washing step, 100 microlitres/well of mouse anti-human Heart FABP (clone 66E2, HyCult Biotechnology BV, Uden, Netherlands), 0.3
20 microgram/ml in PBS-1%BSA, were incubated for 1h at room temperature (R.T) with shaking. After the washing step, 100 microlitres/well of alkaline phosphatase-labelled anti-mouse immunoglobulin (Dako, Denmark), 1.5mg/ml in PBS, were incubated 1h 30min at room temperature with
25 shaking. After the washing step, 50 microlitres/well of alkaline phosphatase substrate, viz. 1.5mg/ml para-nitrophenylphosphate in diethanolamine, was added and the samples were then incubated for 30min. The reaction was stopped with 100 microlitres/well 1M NaOH. Colour
30 development was measured with a microplate reader at a wavelength of 405nm.

"Blank" assays in buffer were also performed.

- 10 -

CK-MB and Troponin-I measurement

AMI was diagnosed by clinical evaluation and Troponin-I and CK-MB measurements. Samples were centrifuged at 1500g for 15min, and stored at -20°C.

5 Serum CK-MB and Troponin-I levels were determined using a fluorescent microparticle enzyme immunoassay (MEIA) with an automated chemical analyser "AxSYM" system (Abbott Laboratories, Abbott Park, IL, USA). The rate of formation of fluorescent products was directly

10 proportional to the amount of Troponin-I in the sample. The detection limit for Troponin-I was 0.3 micrograms/l. CK-MB measurement is proportional to the amount of fluorescent probes and the detection limit was 0.7 micrograms/l.

15 Statistical analysis

H-FABP levels were expressed in optical densitometry (OD) values either as mean plus or minus SD or as median and inter-quartile range. Troponin-I and CK-MB levels were expressed in concentration units (ng/ml). The non-

20 parametric Mann-Whitney U-test and Kruskal-Wallis H-test were used to compare H-FABP, Troponin-I and CK-MB concentrations in plasma between groups. "PRISM" software was used to elaborate box/whisker and scatter plots. The 95% confidence intervals (CI) and Receiver Operating

25 Characteristic (ROC) curves, defined by "Analyse-it" software for Microsoft "EXCEL", were used to assess the discriminatory time point of the indicators. See Murphy, J.M. et al., "Performance of screening and diagnostic tests", Arch. Gen. Psychiatry 44, 550-555 (1987). $P < 0.05$

30 was considered statistically significant.

ResultsClinical characteristics

Patients from the CJD group were given a complete clinical evaluation. CJD was finally diagnosed with the

- 11 -

help of brain immuno-histology after autopsy. Patients from the Control group were admitted to hospital and CJD and AMI were excluded by clinical evaluation.

Patients from the AMI group were admitted to the
5 hospital with confirmed AMI with high Troponin-I levels (>2ng/ml).

Assay results are shown in Table 1 below.

- 12 -

TABLE 1

Assay type	Control Group plasma	AMI Group plasma	Dementia Group CSF	CJD Group plasma	CJD Group CSF
H-FABP median (25-75%) OD, 405 nm	0.25 (0.23-0.27)	2.89 (2.70-3.0)	0.20 (0.16-0.31)	0.79 (0.74-0.86)	0.46 (0.38-0.54)
Troponin-I median (25-75%) IU ng/ml	0 (0.0-0.0)	50 (50-359)	0 (0.0-0.2)	0 (0.0-0.2)	0 (0.0-0.2)

H-FABP plasma levels (OD measurement) in the AMI group were significantly higher than the respective level in the Control group (Table 2). The AMI group had a H-FABP median level (range 25-75%) of 2.89 (2.70-3.0) while the Control group had a level of 0.25 (0.23-0.27). The H-FABP plasma level in the CJD group was between the slopes of the AMI and the Control groups. H-FABP median (range 25-75%) level in the plasma CJD group was 0.79 (0.74-0.86). The sensitivity, specificity, and predictive accuracy of H-FABP levels beyond the cut off value of 0.30 were 100%, 100% and 100% respectively. To confirm differences in H-FABP concentrations between AMI and Control groups, Troponin-I was assayed. In addition, in order to discriminate AMI and CJD, they were also assayed on CJD samples. The Troponin-I concentration was measured in each group. Troponin-I concentration in the AMI group was significantly ($p > 0.01$) higher than in the Control group.

Discussion

The above results indicate that H-FABP is a potential marker for CJD diagnosis. Since H-FABP was

- 13 -

presented as a marker of acute myocardial infarction a few years ago, CJD and AMI had to be discriminated by another AMI biochemical marker such as Troponin-I or CK-MB. After the discrimination of AMI for CJD patient, the serum as well as the CSF H-FABP concentration could be used as a specific marker of CJD.

In the present study, H-FABP assay allowed a sensitivity, a specificity and a predictive accuracy (OD response > 0.30) of 100%. These values were significantly higher than those obtained in another method of pre-mortem detection of CJD, which makes use of the protein 14-3-3, a dimeric phosphoserine-binding protein. This method involves immunoblotting with anti-14-3-3 antibody. The three dementia patients were positive to anti-14-3-3 immunoblotting. The specificity of 14-3-3 is not limited to CJD but includes also Alzheimer's dementia, cerebral complications from head injury and some other forms of dementia.

Acute myocardial infarction is diagnosed with the help of biochemical marker assays such as cardiac Troponin-I, Creatine-Kinase MB, myoglobin and recently H-FABP assay. The H-FABP level for CJD could interfere with AMI and discrimination between AMI and CJD was made with the use of other AMI markers.

EXAMPLE 2

Samples of plasma or CSF were taken from human patients. The disease from which the patients were suffering was in some cases clearly CJD, either sporadic (sp) or new variant (v), as determined by autopsy. In other cases ("not CJD ?"), the patient has been diagnosed as not having CJD, but since some of these patients are still alive, this has not necessarily been confirmed by autopsy. The samples were assayed for CJD by the anti-

- 14 -

14-3-3 method of the prior art and by the present invention.

5 The anti-14-3-3 immunoblot was carried out by running the samples on a 12% SDS-PAGE gel in tris-SDS-glycine buffer. The proteins were thereafter transferred by semi-dry electroblotting at a constant 200 mA for 3 hours, in CAPS buffer, onto a PVDF membrane. The membrane was blocked, incubated with an anti-14-3-3 polyclonal rabbit IgG antibody from Santa Cruz, Inc. (Cat
10 sc 629, Lot L117), washed with buffer and incubated with the second antibody, a goat anti-rabbit immunoglobulin labelled with horseradish peroxidase (Dako, Denmark). The membrane was then washed again. The washing after each incubation was done in PBS buffer, pH 7.2, with 5%
15 "Tween" three times quickly and five times for five minutes each time. The peroxidase was then assayed by a standard enhanced chemiluminescence method, using a Boehringer Mannheim kit, "BM Chemiluminescence Blotting Substrate (POD)". The luminescence observed denoted a
20 positive result in the immunoblotting.

The method of the present invention was as described in Example 1, except that the sensitivity cut-off applied (using ROC curves) was at OD >0.2 for plasma samples and OD >0.1 for CSF samples. Table 2 shows the results.

25 Referring to Table 2, the anti-14-3-3 test was performed twice, by different operatives in the inventors' laboratory, yielding the same results. The correlation between the anti-14-3-3 and the anti-H-FABP results was nearly 100%, the exception being the sample
30 CSF-10, where the result was not clear. The plasma samples gave positives with anti-H-FABP in four cases in which the anti-14-3-3 test gave a negative. This could mean that the anti-14-3-3 test is not giving a true result in all cases.

35

TABLE 2

Sample Designation	Disease Assignment	Anti-14-3-3 Immunoblot (Prior art)*	Anti-H-FABP ELISA (This inv.)
PLAS2	vCJD	Negative	Positive
PLAS3	vCJD	Negative	Negative
PLAS4	vCJD	Negative	Positive
PLAS5	spCJD	Positive	Positive
PLAS6	spCJD	Negative	Negative
PLAS7	spCJD	Positive	Positive
PLAS9	not CJD ?	Positive	Positive
PLAS10	not CJD ?	Positive	Positive
PLAS11	not CJD ?	Negative	Positive
PLAS12	not CJD ?	Negative	Positive
CSF1	spCJD	Positive	Positive
CSF2	spCJD	Positive	Positive
CSF3	spCJD	Positive	Positive
CSF4	spCJD	Positive	Positive
CSF5	spCJD	Positive	Positive
CSF10	vCJD	Positive	Positive
CSF11	vCJD	Positive	Unclear
CSF12	vCJD	Positive	Positive
CSF6	not CJD ?	Negative	Negative
CSF7	not CJD ?	Positive	Positive
CSF8	not CJD ?	Negative	Negative
CSF9	not CJD ?	Negative	Negative
CSF13	not CJD ?	Negative	Negative
CSF14	not CJD ?	Negative	Negative

- 16 -

* Performed twice, by different workers, with the same results.

- 17 -

EXAMPLE 3

The method of the invention was carried out on pooled, concentrated, samples of CSF from 4 cattle diagnosed as having BSE and on pooled, concentrated
5 samples from 3 healthy cattle as controls. (The samples were concentrated with "Microcon" ,from Amicon, in order to increase the signal to background ratio).

A rat/mouse H-FABP ELISA kit from Hycult Biotechnology B.V., Uden, The Netherlands, was used,
10 according to the manufacturer's instructions, the assay being similar in principle to the sandwich ELISA described in Example 1. However, the first antibody, bound to the wells, was anti-rat/mouse H-FABP, rather than anti-human H-FABP, and the second antibody was
15 peroxidase-labelled, anti-rat/mouse. (These antibodies appear to be anti- to both rat and mouse. It should be explained that this kit was not intended to detect bovine H-FABP. It was found unexpectedly in the present invention that the anti-rat/mouse H-FABP antibody
20 recognises bovine H-FABP). The assay is colorimetric, using SMP substrate and with readout at 450 nm.

The results, shown in Table 3, are the average of duplicate assays and indicate clearly the difference
25 observed in the BSE-affected cattle compared with the healthy cattle.

- 18 -

TABLE 3

SAMPLE	Average intensity	Coefficient of variation
Blank (PBS)	0.172	3.6 %
Healthy CSF	0.178	11.8 %
Healthy CSF	0.189	2.4 %
BSE CSF	0.304	1.5%
BSE CSF	0.576	4.0%
BSE CSF	0.465	10.8%
Bovine heart (10 mg/ml.)	2.872	2.0%
Blank (PBS)	0.178	2.1%

* * * * *

- 5 Each of the above cited publications is herein incorporated by reference to the extent to which it is relied on herein.

- 19 -

CLAIMS

1. A method of diagnostic assay for a transmissible spongiform encephalopathy (TSE) or the possibility thereof in a sample of body fluid taken from a subject suspected of suffering from the TSE, which comprises determining the concentration of heart or brain fatty acid binding protein (H-FABP or B-FABP) in the sample.
2. A method according to Claim 1, wherein the subject is a human patient and the concentration of H-FABP is determined in a first assay, whereby a positive result indicates either a CJD or acute myocardial infarction, and which further comprises carrying out a second diagnostic assay, for acute myocardial infarction (AMI) only, whereby a positive result in the H-FABP assay and a negative result in the assay for AMI indicates that the patient is or might be suffering from a CJD.
3. A method according to Claim 2, wherein the assay for AMI comprises determining the concentration of troponin-1 or creatine kinase MB in plasma.
4. A method according to Claim 1, 2 or 3, wherein an antibody to H-FABP is used in the assay for H-FABP.
5. A method according to Claim 4, wherein the subject is a human patient and a mouse anti-human FABP monoclonal antibody is used.
6. A method according to Claim 4 or 5, wherein the assay for H-FABP comprises a sandwich ELISA.
7. A method according to Claim 1, wherein B-FABP or an antibody thereto is used without any assay for AMI in combination therewith.
8. A method according to any preceding Claim, wherein the H-FABP or B-FABP assay is carried out on a blood plasma or serum sample.

1/1




Fig. 1